# *testoval*<sup>®</sup>

## CHROMAT-TEST

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### Directions for Use

Remove the test cell from the comparator and fill it up to the mark with water sample (turbid samples must be cleared by filtration). Add 4 drops of reagent A and 8 drops reagent B, mix by stirring with the red spoon.

Reinsert the test cell into the comparator and after 3 minutes hold the comparator against a bright background or light source, and match the colour intensity of the central field (test solution) with one of the colours of the peripheral colour fields. Read the value in ppm Cr, printed below the matching colour field.

Should the colour in the central field be more intensive than any of the colours in the peripheral fields, repeat the test with a water sample diluted 1:10. To effect the dilution, fill the graduated plastic tube up to the lower mark with water sample, then add chromate-free water (e. g. distilled water) up to the upper mark. Transfer this solution into the test cell which previously has been rinsed with chromate-free water and carry out the test as described above. The value obtained must in this case be multiplied by 10.

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