testoval[®]

PHOSPHATEST[®]

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1 ppm P₂O₅ \triangleq 1,34 g/m³ PO₄³⁻

Directions for Use

The molybdenum-complex method, on which practically all colorimetric phosphate determinations are based, has been significantly improved for this phosphate test by the addition of other complex forming species and a novel reducing agent. Therefore the colour development is complete only 5 minutes after addition of the reagents to the sample. This compares favourably with other similar methods where colour changes still occur even after 10 minutes and where sometimes it takes 30 minutes until the colour is sufficiently stable to be evaluated.

Remove the test cell from the comparator and fill it up to the mark with water sample (turbid samples must be cleared by filtration). Add 4 drops each of reagent A and reagent B and 8 drops reagent C, mix by stirring with the red spoon.

Reinsert the test cell into the comparator and after 5 minutes hold the comparator against a bright background or light source, and match the colour intensity of the central field (test solution) with one of the colours of the peripheral colour fields. Read the value in ppm P₂O₅, printed below the matching colour field.

Should the colour in the central field be more intensive than any of the colours in the peripheral fields, repeat the test with a water sample diluted 1:10. To effect the dilution, fill the graduated plastic tube up to the lower mark with water sample, then add phosphate-free water (e. g. distilled water) up to the upper mark. Transfer this solution into the test cell which previously has been rinsed with phosphatefree water and carry out the test as described above. The value obtained must in this case be multiplied by 10.

Should the colour in the central field even in the 1:10 diluted sample be more intensive than any of the colours in the peripheral fields, the phosphate contents of the sample is excessively high, i. e. over 100 ppm.

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