SILICAT-TEST

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SILICATE TEST

Directions for Use

Remove the test cell from the comparator and fill it up to the mark with the water sample.

Add 2 drops of reagent A and 2 drops of reagent B. Stir with the red spoon.

After 3 minutes, add 5 drops of reagent C and stir. After a further minute, add a pipette of reagent D filled to the mark and stir once again. Insert the test cell into the comparator. After 15 minutes, hold the comparator against a light source or bright background and match the colour of the central field (test solution) with one of the colours of the peripheral colour fields. Read the value shown under each matching colour field.

Should the colour in the central field be more intensive than the highest colour value in the peripheral fields, dilute the water sample 1:10. For this purpose, fill the graduated plastic tube up to the lower mark with the water sample and dilute with silicate-free water (e.g. distilled water) to the upper mark. Transfer the diluted sample into the test cell which has been rinsed previously with silicate-free water and proceed as described above. The value obtained is then multiplied by 10.

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