


pH Measurement To Determine Acidification of Sushi Rice

The rice used for sushi must be acidified with acetic acid (vinegar) to pH less than 4.6 to inhibit the growth of pathogenic bacteria. To measure pH, simply place a sample of rice mixture onto the flat sensor of LAQUAtwin pH meter.



Scan QR code with your mobile device to access online content



LAQUAtwin pH Meters Product Page

Introduction

Rice used for sushi must have a pH of less than 4.6. At pH levels below 4.6, most pathogenic bacteria do not grow or produce toxins¹. Thus, the rice must be acidified using acetic acid (vinegar) to be classified as non-hazardous. The LAQUAtwin pocket pH Meter is used as quality control check to ensure that the rice is adequately acidified, before selling to consumers. This is an easy, quick check method used to abide to the ANZ Standards² in ensuring that customers are safely consuming sushi.

Method

Acetic acid (vinegar) should be mixed into the rice according to the following table.

Ingredients	Recipe 1	Recipe 2	Recipe 3
Short Grain Rice	900 g	900 g	900 g
Water	1100 ml	1320 ml	1250 ml
Rice Vinegar	135 ml	99 ml	128 ml
Sugar	57 g	94 g	44 g
Salt	9 g	25 g	8 g

A small sample of the rice mixture is placed on the flat sensor of the LAQUAtwin pocket pH Meter and measured. If the measured value is above pH 4.6, add more acetic acid to the rice mixture and stir well. Place new rice sample on the sensor and repeat testing process. After tests, wash the sensor with diluted soap water and pat dry with a paper tissue.

Results and Benefits

The use of accurate pH testing in controlling the quality of sushi rice prevents the growth of pathogenic bacteria and toxins.


The LAQUAtwin pocket pH meter is small and compact; convenient to carry around in your pocket and is ideal for on-site testing. Its easy-to-use interface makes the LAQUAtwin pocket pH Meter an indispensable tool for food testing. 💧

¹Hocking, A.D; 2003. Foodborne Microorganisms of Public Health Significance, AIFST, Waterloo
²Food Safety Guideline for Preparation and Display of Sushi, June 2007, NSW/FA/F1005/0706

LAQUAtwin

Pocket ION Meter

Unique Features



LAQUAtwin: the only meters with flat sensor technology.

HORIBA's highly-sensitive, flat sensor technology opens up new possibilities for sampling and sample types. Only a small amount of sample is required, so you can easily sample in situ without the need for beakers or other labware. Sensors are easily replaced as required.

Calibrate and measure at the touch of a button—the smiley face will tell you when the result can be read.

Hassle-free automatic calibration with a few drops of standard solution reassures you of your measurement accuracy. Two-point calibration is also possible.*1

*1 Except for B-711



LAQUAtwin is fully waterproof and dustproof.

The meter and sensor are fully waterproof*3 and dustproof, so you can take it anywhere.

*3 IP67 rated. Will withstand immersion for 30 minutes at 1 m. Not suitable for underwater use.

Carry case comes as standard for handy portability.


The compact carry case contains everything you need for your measurements, including the standard solution and sampling sheets.



1 X 6


One meter, six methods.

Only LAQUAtwin allows you to be this flexible!
Choose the best method according to your sample, your situation, and your needs.



01 Immersion

When you're in the lab, you can test the sample in a beaker. Ensure the sensor guard sliding cap is open.




02 Scoop

Use as a scoop to test water eg from a river. A vertical scoop for an aquarium is also available with a unique sensor guard.



03 Drops

Place a drop of the sample onto the sensor with a pipette. LAQUAtwin meters can measure sample volume as low as 0.1mL



04 Solid Samples

Foods containing some moisture can be tested by placing a small piece directly onto the sensor.



05 Powders

LAQUAtwin meters can also test dry powders. Simply place the powder sample onto the sensor and drop on your defined volume of pure water.



06 Paper and textiles

To test sheets of paper and textiles, cut up the sample into small pieces and place directly onto the sensor. Drop on your defined volume of pure water.

LAQUAtwin Pocket Ion Meters Lineup



HORIBA Instruments (Singapore) Pte. Ltd.
 83 Science Park Drive, #02-02A, The Curie,
 Singapore 118258
 Phone: 65 6908-9660
 Fax: 65 6745-8155
www.horiba-laqua.com
 e-mail: laqua@horiba.com

HORIBA UK Limited
 Kyoto Close, Moulton Park, Northampton
 NN3 6FL
 Phone: 44 (0) 1604 542567
 Fax: 44 (0) 1604 542699
www.horiba.com/uk
 e-mail: waterquality@horiba.com

HORIBA Instruments Incorporated
 9755 Research Drive, Irvine, California
 92618 USA
 Phone: +1 949 250 4811
 FAX: +1 949 250 0924, +1 949 468 1890
www.horiba.com/us/en



www.horiba-laqua.com

IMS
 HORIBA Group is operating Integrated
 Management System (IMS)
 ISO9001 JQA-0298 / ISO14001
 JQA-E-90039 / ISO13485
 JQA-MD0010 / OHSAS18001 JQA-
 OH0068Rorum perid novis esimoente,



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA