
Product	9618S-10D
Part number	3200585447
pH Range	0 to 14
Temperature Range	0 to 60
Liquid Junction	Ceramic
Internal electrode	Silver/silver chloride
Reference electrode internal solution	3.33 mol/L-KCl solution (model 300)
Temperature Sensor	Integrated
Wetted material	Glass, ceramic
Cable	1m
Connectors	BNC, phono jack
Overall Length	151.5mm
Diameter	3mm, 12mm
Filling Solution	525-3 or 300 (3.33M KCl)
Waterproof	Yes

Description

This glass-body, refillable pH electrode with temperature sensor can measure samples as low as 50µL.

Compatible with extremely small sample containers (e.g., microtubes)

Temperature sensor is placed beside the bulb for quick response



Dimensions in mm

